

Implants with conical connection

Product overview



Achieve high function and esthetics with the internal conical connection



Access innovative restorative solutions


High predictability with the cement-free angulated screw channel abutment, or the full-contour zirconia implant crown. Enjoy easy handling with the unique Omnigrip tooling's pick-up function and its incredible grip on the screw.

Create natural-looking esthetics

Built-in platform shifting is designed to increase soft tissue volume.

Benefit from a strong sealed connection

Internal conical connection with hexagonal interlocking provides a tight seal and high mechanical strength.



One restorative connection, three implant designs

Using one restorative connection, you benefit from the same restorative assortment and routine, while having the surgical flexibility of three implant designs that each provide high initial stability.



NobelParallel Conical Connection – the parallel-walled design

- Proven and documented parallel-walled implant body for universal use in all bone qualities and a wide range of indications.
- Straightforward and flexible surgical protocol shortens treatment times.
- High initial stability provided by the implant design and surgical protocol enables Immediate Function.



NobelReplace Conical Connection – the tapered design

- Proven tapered implant body for high initial stability with all loading protocols.
- Efficient treatment flow with step-by-step tapered drilling protocol and consistent color-coding of implants, drills and prosthetic components.



NobelActive – the hybrid design

- Unique combination of straight drilling protocol and tapered implant design provides exceptionally high initial stability, even in compromised bone situations.
- Back-tapered collar for stable bone levels and high esthetics.

Order your implant of choice today

NobelParallel Conical Connection*

Platform	Implant Ø	Length						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
NP	3.75 mm	37963	37964	37965	37966	37967	37968	37969
RP	4.3 mm	37970	37971	37972	37973	37974	37975	37976
RP	5.0 mm	37977	37978	37979	37980	37981	37982	37983
WP	5.5 mm	37984	37985	37986	37987	37988	37989	-



NobelActive

Platform	Implant Ø	Length						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
3.0	3.0 mm	-	-	36769	36770	36771	36772	-
NP	3.5 mm	-	35221	34125	34126	34127	34128	35215
RP	4.3 mm	-	35223	34131	34132	34133	34134	35219
RP	5.0 mm	-	35225	34137	34138	34139	34140	35220
WP	5.5 mm*	37806	37807	37808	37809	37810	37811	-



NobelReplace Conical Connection

Platform	Implant Ø	Length				
		8 mm	10 mm	11.5 mm	13 mm	16 mm
NP	3.5 mm	36699	36700	36701	36702	36703
RP	4.3 mm	36704	36705	36707	36708	36709
RP	5.0 mm	36710	36711	36712	36713	36714



With TiUnite on collar


NobelReplace Conical Connection PMC (Partially Machined Collar)*

Platform	Implant Ø	Length				
		8 mm	10 mm	11.5 mm	13 mm	16 mm
NP	3.5 mm	37284	37285	37287	37288	37289
RP	4.3 mm	37290	37291	37292	37293	37294
RP	5.0 mm	37295	37296	37297	37298	37299



With 0.75 mm machined collar

Cover screw

	3.0	NP	RP	WP
	36775	36649	36650	37812

* Cover screw included.



On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit nobelbiocare.com/warranty



The key to successful treatments

Diagnose and plan implants in 3D based on the patient's anatomy and prosthetic needs with the NobelClinician Software. For more complex cases, you have the option to place the implants fully guided using the well-established NobelGuide concept.

nobelbiocare.com/nobelclinician

Surgical components

Shared components

Precision Drill		Guide Drill	
3.0, NP, RP, WP	36118	NP, RP, WP	35426




Implant Drivers	28 mm	37 mm
3.0	36773	36774
NP	36718	36719
RP	36720	36721
WP	37859	37860




Bone Mills and Guides	Bone Mill with Guide	Bone Mill Guide
3.0	37861 Ø4.0	37862
NP	37863 Ø4.4, 37864 Ø5.2	37865
RP	37866 Ø5.2, 37867 Ø6.2	37868
WP	37869 Ø6.7	37870



37888 Bone Mill Kit Conical Connection
Includes Bone Mills with Guide for all conical connection platforms.



NobelActive 3.0 for limited anterior spaces.




NobelActive WP and NobelParallel CC WP with NobelProcera Full-Contour Zirconia implant crown (0–25° angulation) is ideal for cases in the posterior region.



Solutions for all indications

For the full assortment of Nobel Biocare products, see our current product catalog.

nobelbiocare.com/catalog



Order online

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

nobelbiocare.com/store



Order by phone or fax

Call our customer service team, send a fax or contact your sales representative.

Straight drilling protocol

NobelActive and NobelParallel Conical Connection

Twist Drills	7–10 mm	7–15 mm	10–18 mm
Ø 1.5	–	31278	–
Ø 2	32296	32297	32299



Twist Step Drills	7–10 mm	7–15 mm	10–18 mm
Ø 2.4/2.8	32260	32261	32262
Ø 2.8/3.2	37873	34638	34639
Ø 3.2/3.6	32363	32264	32265
Ø 3.8/4.2	32275	32276	32277
Ø 4.2/4.6	37874	34582	34583
Ø 4.2/5.0	37875	37876	–



Cortical Drills NobelParallel CC		
Ø 3.75		38000
Ø 4.3		38001
Ø 5.0		38002
Ø 5.5		38003



Screw Taps NobelParallel CC	7–10 mm	7–13 mm	7–15 mm	7–18 mm
Ø 3.75	–	37990	–	37991
Ø 4.3	–	37992	–	37993
Ø 5.0	–	37994	–	37995
Ø 5.5	37996	–	37997	–



Tapered drilling protocol

NobelReplace Conical Connection (and PMC)

Drill with Tip Tapered Ø 2 mm	
NP, RP	36117




Drills	8 mm	10 mm	11.5 mm	13 mm	16 mm
NP 3.5 mm	32075	29367	36113	29368	29369
RP 4.3 mm	32076	29370	36114	29371	29372
RP 5.0 mm	32077	29373	36115	29374	29375



Dense Bone Drills	13 mm	16 mm
NP 3.5 mm	29377	29378
RP 4.3 mm	29380	29381
RP 5.0 mm	29383	29384




Screw Taps	
NP 3.5 mm	36717
RP 4.3 mm	32090
RP 5.0 mm	32091



37883 NobelActive Surgery Kit

Includes instruments to perform implant surgery with NobelActive implants. Drills and screw taps must be ordered separately.



38004 NobelParallel CC Surgery Kit

Includes instruments to perform implant surgery with NobelParallel Conical Connection implants. Drills and screw taps must be ordered separately.



37465 NobelReplace CC Surgery Kit

Includes instruments and drills to perform implant surgery with all NobelReplace Conical Connection implants.





36915 Implant Driver Kit Conical Connection



Complementary kit for 36894 NobelReplace Tapered Surgery Kit. Includes implant drivers for NP and RP implants with conical connection.







Prosthetic components





Healing abutments

	3.0			NP			RP			WP		
	Healing Abutment	3 mm	5 mm	7 mm	3 mm	5 mm	7 mm	3 mm	5 mm	7 mm	3 mm	5 mm
	Ø3.2	36794	36795	36796	-	-	-	-	-	-	-	-
	Ø3.6	-	-	-	36639	36640	36867	36643	36644	36872	-	-
	Ø3.8	36797	36798	36799	-	-	-	-	-	-	-	-
	Ø5.0	-	-	-	36641	36642	36868	36645	36646	36873	37813	37814
	Ø6.0	-	-	-	-	-	-	36647	36648	36874	-	-
	Ø6.5	-	-	-	-	-	-	-	-	-	37815	37816
	Slim Healing Abutment	-	37669	37670	-	37666	37665	-	37667	37668	-	-

		NP			RP			WP	
	Healing Abutment Bridge	3 mm	5 mm	7 mm	3 mm	5 mm	7 mm	3 mm	5 mm
	Ø4.0	36864	36865	36866	-	-	-	-	-
	Ø5.0	-	-	-	36869	36870	36871	-	-
	Ø6.0	-	-	-	-	-	-	37817	37818
	Healing Abutment Anatomical PEEK				WP				
					6x7 mm		7x8 mm		
					37819		37820		


Temporary abutments¹


		3.0		NP		RP		WP	
	Immediate Temporary Abutment²	1.5 mm	3 mm	1.5 mm	3 mm	1.5 mm	3 mm	1.5 mm	3 mm
		36777	36778	36653	36655	36654	36656	-	-
	Temporary Abutment Engaging	36779	-	36663	-	36664	-	37823	37824
	Temporary Coping Multi-unit Titanium⁴	-	-	29046	-	29046	-	29046	-
	Slim Temporary Abutment	6.5 mm	7.5 mm	6.5 mm	7.5 mm	6.5 mm	7.5 mm	-	-
		37675	37676	37671	37672	37673	37674	-	-


		NP		RP		WP	
	QuickTempTM Abutment³	1.5 mm	3 mm	1.5 mm	3 mm	1.5 mm	3 mm
		36659	36657	36660	36658	-	-
	Temporary Abutment Non-Engaging	36661	-	36662	-	37825	37826
	Temporary Coping Multi-unit Plastic⁵	DCA	468-0	DCA	468-0	DCA	468-0
	Temporary Abutment Anatomical PEEK	-	-	-	-	6x7 mm	7x8 mm
		-	-	-	-	37821	37822

Impression copings




Implant level




		3.0	NP	RP	WP	
	Closed Tray	13 mm	13 mm	9 mm	13 mm	9 mm 13 mm
	Ø3.3	36801	-	-	-	-
	Ø3.6	-	36538	36541	36540	-
	Ø3.8	36803	-	-	-	-
	Ø5.0	-	36539	36543	36542	37851 37850
	Ø6.0	-	-	36545	36544	-
	Ø6.5	-	-	-	-	37853 37852

		3.0	NP	RP	WP	
	Open Tray	14 mm	10 mm 14 mm	10 mm 14 mm	10 mm 14 mm	10 mm 14 mm
	Ø3.2	36800	-	-	-	-
	Ø3.6	-	36258	36260	36263	36262
	Ø3.8	36802	-	-	-	-
	Ø5.0	-	36259	36261	36265	36264
	Ø6.0	-	-	-	36267	36266
	Ø6.5	-	-	-	-	37857 37856

	Bridge Open Tray	-	36930	36931	37858
---	-------------------------	---	-------	-------	-------


Abutment level

		NP	RP	WP
	Closed Tray Multi-unit	29090	29090	29090
	Bar Closed Tray Multi-unit	29093	29093	29093
	SnappyTM Abutment			
	4.0	35995	35996 35997 Wide	37841
	5.5	35868	35869 35870 Wide	37846

		NP	RP	WP
	Open Tray Multi-unit⁶	29089	29089	29089
	Locator[®] Abutment	4/pkg REF08505		20/pkg REF08505-20
	Block Out Spacer	-	-	REF08514

Prefabricated final abutments¹

Zirconia


	NP	RP
	Procera® NP #1 36627	RP #5 36631
	NP #2 36628	RP #6 36632
	Esthetic Abutment ⁷ NP #3 36629	RP #7 36633
	NP #4 36630	RP #8 36634
	RP #9 36635	
	RP #10 36636	
	RP #11 36637	
	RP #12 36638	


Consistent color-coding


Prosthetic components are color-coded for accurate and fast identification.





Titanium


	3.0		NP		RP		WP	
	Esthetic Abutment	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	6x7 mm		
		36782 36783	36665 36666	36669 36671	37827			
		4.5 mm 36814	4.5 mm 36249	4.5 mm 36251	7x8 mm 37828			


	Narrow Profile Abutment	7 mm 9 mm	7 mm 9 mm	7 mm 9 mm	7 mm 9 mm			
		36781 36780	36678 36679	36680 36681				

	NP		RP		WP	
	Snappy™ Abutment 4.0 ⁸	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	
		36695 36696	36693 36694	37834 37835		
			36691 Wide 36692 Wide			

	Snappy™ Abutment 5.5 ⁸	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	
		36684 36688	36682 36683	37836 37837		
			36689 Wide 36690 Wide			


	3.0		NP		RP	
	Esthetic Abutment 15°	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	
		36784 36785	36667 36668	36672 36673		
		4.5 mm 36815	4.5 mm 36250	4.5 mm 36252		


	Locator® Abutment					
	1.0 mm	–	REF02069	REF02075		
	2.0 mm	–	REF02070	REF02076		
	3.0 mm	–	REF02071	REF02077		
	4.0 mm	–	REF02072	REF02078		
	5.0 mm	–	REF02073	REF02079		
	6.0 mm	–	REF02074	REF02080		

	NP		RP		WP	
	Universal Base Abutment	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	1.5 mm 3 mm	
		38213 38216	38214 38217	38215 38218		


	Universal Base burn-out coping	38221	38221	38222		
---	--------------------------------	-------	-------	-------	--	--


Gold


	3.0	NP	RP	WP	
	GoldAdapt™ Engaging	36790	36728	36729	37849


	NP	RP	WP	
	GoldAdapt™ Non-Engaging	36726	36727	37848


Multi-unit


	NP	RP	WP	
	Multi-unit Abutment			
	1.5 mm	36611	36616	37829
	2.5 mm	36613	36617	37830
	3.5 mm	36624	36625	37831
	4.5 mm	–	36626	–

	Healing Cap Multi-unit			
	1/pkg	31145	31145	31145
	5/pkg	29064	29064	29064
	Wide 1/pkg	31146	31146	31146

	NP		RP		WP	
	17° Multi-unit Abutment	2.5 mm 3.5 mm	2.5 mm 3.5 mm	2.5 mm 3.5 mm	2.5 mm 3.5 mm	
		36614 36615	36618 36619	37832 37833		

	30° Multi-unit Abutment	3.5 mm 4.5 mm	3.5 mm 4.5 mm		
		36620 36621	36622 36623		

	Gold Coping Multi-unit ⁴			
	1/pkg	29043	29043	29043
	5/pkg	29042	29042	29042

	Gold Coping Bar Multi-unit ⁴	29045	29045	29045
---	---	-------	-------	-------

1 Abutment screw included with all abutments.
 2 Plastic coping 31656 included (cannot be used as a burn-out).
 3 Plastic coping 33404 included (cannot be used as a burn-out).
 4 Prosthetic screw included.
 5 Prosthetic screw not included (cannot be used as a burn-out).
 6 15 mm guide pin included.

7 Prefabricated zirconia abutments with internal conical connection are not indicated for posterior use.
 8 Includes abutment screw, impression coping, healing cap and plastic/temporary coping engaging (plastic/temporary coping non-engaging also available).





NobelProcera CAD/CAM for restorative flexibility and predictable results



- Titanium or zirconia abutments that support the surrounding tissue with an individualized emergence profile
- Cement-free angulated screw channel (ASC) zirconia abutments for better esthetics, optimal occlusal function and easier access
- Cement-free and predictable full contour zirconia implant crown with angulated screw channel
- Implant bridges in zirconia and titanium, as well as fixed and fixed-removable implant bars overdenture for improved quality of life

Contact your local Nobel Biocare sales office to learn more about the full range of precision-milled NobelProcera restorations.











Replicas

	3.0	NP	RP	WP
 Implant	36791	36697	36698	37879
	NP	RP	WP	
 Snappy™ Abutment				
4.0	35982	35983	35984 Wide	37842
5.5	35854	35855	35856 Wide	37847

	NP		RP		WP	
 Multi-unit Abutment	1/pkg	5/pkg	1/pkg	5/pkg	1/pkg	5/pkg
	31161	29110	31161	29110	31161	29110
	4/pkg		20/pkg			
 Locator® Abutment	Ø4 mm		REF08530		REF08530-20	
	Ø5 mm		REF08516		REF08516-20	











Screws









	3.0	NP	RP/WP
 Clinical Screw	37890	37891	37892
 Omnigrip Clinical Screw*	-	37367	37606
 Screw Multi-unit Angled	-	36892	37893
 Prosthetic Screw Multi-unit	-	29285	29285

	3.0	NP	RP/WP
 Laboratory Screw	36805	37894	37895 (5/pkg)
 Omnigrip Laboratory Screw*	-	37374	37607
 Laboratory Screw Multi-unit Angled	-	37896	37897
 Laboratory Screw Multi-unit 5/pkg	-	29287	29287

* for NobelProcera ASC abutments and FCZ Implant crowns

Other components

 37448 Prosthetic Kit			
Unigrip™ Screwdriver Machine 20 mm	29151		
Unigrip™ Screwdriver Machine 30 mm	29153		
Multi-unit Screwdriver Machine 21 mm	29158		
Omnigrip Screwdriver Machine 20 mm	37379		
Omnigrip Screwdriver Machine 30mm	37381		
Manual Torque Wrench Prosthetic	29165		
Implant Driver for Slim Abutment			
Implant Driver CC 3.0 for Slim Abutment	37713		
Implant Driver CC NP for Slim Abutment	37677		
Implant Driver CC RP for Slim Abutment	37678		

 37508 Abutment Retrieval Kit			
Abutment Retrieval Instrument Zirconia CC NP	37512		
Abutment Retrieval Instrument Zirconia CC RP/WP	37882		
Abutment Retrieval Tool Titanium NP	36247		
Abutment Retrieval Tool Titanium RP/WP	37881		
Abutment Release Pin 3.0	37509		
Abutment Release Pin NP	37510		
Abutment Release Pin RP/WP	37511		



nobelbiocare.com